



ADF- B180 Blood Pressure meter

Application note for App developer

Device Name : BPM

Service UDID: FFF0

Characteristics UDID : FFF4

Measurement Data format:

Byte No.	Type	Description	Unit
1	Flag	● Bit0 (Unit): ■ 0 = mmHg ■ 1=KPA ● Bit1 (Time Stamp): ■ Default = 0 ● Bit2 (Pulse Rate) ■ 0 = Without ■ 1=With Pulse Rate ● Bit3 (User ID) ■ 0 = Without ID ■ 1 = With User ID ● Bit4 (Meas. Status) ■ 0 = Without status ■ 1 = With status ● Bit5 (Test Status) ● 0 = Normal data ● 1 = Testing data (if bit 5 = 1 and other bit = 0. Please see Testing data format)	
2	Systolic – High	Systolic data	mmHg
3	Systolic – Low		
4	Diastolic – High	Diastolic data	mmHg
5	Diastolic – Low		

6	Map – High	Default = 0x00	
7	Map – Low		
8	Pulse Rate – High	Pulse rate data	Pulse per min.
9	Pulse Rate – Low		
10	User ID	01, 02	
11	PAD	1=PAD detected	
12	Meas. Status	0=normal 01=E-1 Pulse and pressure measurement fail 02=E-2 Noise Interference 03=E-3 Inflation fail 05=E-5 abnormal BP 0B=E-B Battery Low	

Testing Data(current measurement pressure status) format:

Byte No.	Type	Description	Unit
1	Flag	Default = 0x20	
2	Pressure	Current Pressure data	mmHg

Transfer Memory format:

Byte No.	Type	Description	Unit
1,2	Header	FD FD	
3	Status	Last memory : FF No memory: FF Etc: FE	
4	User	1~2	
5~7	Data	SYS 0~255 DIA 0~255 PUL 0~255 if no memory all 0	
8~13	Date and Time	Year1 20 Year2 19~99	

		Month 1~12 Day 1~31 Hour 0~23 Minute 0~59 if no memory all 0	
14~15	Footer	0D 0A	

Service UDID: FFF0

Characteristics UDID : FFF3

Update Time format:

Byte No.	Type	Description	Unit
1-4	Header	FD FD FA 09	
5-10	Date and Time	Year 19~99 Month 1~12 Day 1~31 Hour 0~23 Minute 0~59 Second 0~59	
11~12	Footer	0D 0A	

Start Measurement:

Byte No.	Type	Description	Unit
1-4	Header	FD FD FA 05	
5~6	Footer	0D 0A	

Power Off:

Byte No.	Type	Description	Unit
1-4	Header	FD FD FA 06	
5~6	Footer	0D 0A	



Transfer Memory:

Byte No.	Type	Description	Unit
1-4	Header	FD FD FA 07	
5~6	Footer	0D 0A	

Received Memory ACK:

Byte No.	Type	Description	Unit
1-4	Header	FD FD FA 60	
5~6	Footer	0D 0A	

Remark: APP send received memory command to BLE and device will clear all device memories.

